

Tips and Tricks to make your pumping journey with Magic InBra™ simple

1. Start by reading the instructions for use and watching the how-to use videos. Found here [Medela Magic InBra™](#)
2. Make sure you assemble all parts correctly. Focus on the membrane, valves and container lid assembly. The more you use your pump, the easier this gets.
3. If the pump stops early into the pumping session:
 - a. Check the assembly of all parts. When disassembling make sure the pump is first completely turned off.
 - b. If the issue persists, try using a smaller shield size to engage the breast into the pump more effectively. We have 21mm and 24mm shields included in the box. 18mm and 27mm are available for purchase. If you feel unsure about purchasing a new size, please contact our customer service for guidance.
4. If milk leaks, make sure the container lid is closed firmly and securely. Also make sure you pause the pump and drain the milk into the container when the milk fill level reaches the 'max' line/130 ml. Remember, around 20 ml still needs to be drained from the tunnel into the container.
5. We urge you to always pump at a comfortable vacuum level. If pumping feels painful, follow the instructions to find your maximum comfort vacuum. The pump will remember this setting for subsequent pumping sessions. It is also important to reassess this vacuum level every few weeks or when needed.
6. If you are still having trouble using the pump and the app. We are here to help! **Scan the QR code below and click on 'Contact us' to contact your local customer service for personalized assistance.**

