Thopaz+™

Cardiothoracic Drainage

EN Instructions for use

DA Brugsanvisning

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NO Bruksanvisning

SV Bruksanvisning

HR Upute za upotrebu

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SR Uputstva za upotrebu

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Table of Contents

Congratulations	6
1 Warnings and safety instructions	7
Safety instructions	11
2 Power supply and battery operation	12
Backup battery	12
Charge battery	12
Battery low	13
Battery empty	13
3 Description	14
Introduction	14
Intended use/purpose	14
Indications for use	14
Contraindications	14
Intended user	14
Intended patient population	14
Important note	14
Undesirable side effects	14
4 Overview	15
Main elements of Thopaz+	15
Display	16
How to activate a button	16
Tubings	17
Canisters	19
5 Additional information	20
Definition of vacuum	20
Regulated pressure	20
Safety pressure limit	20
Tube flushing	20
Safety chamber in canister	20
Fluid measurement	21
Limitations	21
Overflow protection/bacteria filter	22
Light sensor	22
6 Installation	23
Check initial delivery	23
Initial startup	24
7 Preparation for use/operating instructions	26
Checks before use	26
Connect tubing	26
Snap in canister	27
Using 2l canister	27
Switch Thopaz+ on	28

Carry out functional check.....	29
Check therapy progress.....	30
Main display.....	30
Reset fluid display.....	31
History.....	32
Catheter check.....	34
Switching Thopaz+ off.....	35
Transferring data to PC with ThopEasy+.....	36
8 Replace canister.....	37
Canister with solidifier: activation of the solidification process.....	38
9 Changing pressure during operation.....	39
Physiological pressure.....	39
10 Changing settings.....	40
Factory settings.....	40
Change settings.....	41
Save a changed parameter as the new standard.....	41
11 Taking a drainage sample.....	42
12 Troubleshooting and information signals.....	43
Description.....	43
13 General reprocessing guidelines.....	46
14 Warranty, maintenance and checks.....	50
Warranty.....	50
Maintenance.....	50
Routine check.....	50
Instructions for FPT test.....	50
Safety check.....	51
15 Disposal.....	52
16 Accessories.....	52
17 Technical specifications.....	53
18 Meaning of symbols.....	54
19 Technical documentation (English only).....	472
Appendix A.....	478
Appendix B.....	480

These instructions for use belong to firmware 1.01.

Congratulations

With Thopaz+, you have acquired an innovative digital chest drainage system that sets new standards in therapy and chest drainage management.

Thopaz+ has an electronic measuring and monitoring system with optical and acoustic status indications. The device is a dry system, which means no fluids are necessary for operation. Important information regarding the course of therapy is indicated digitally and graphically on the display. Upon completion of therapy, it can be transferred to a PC.

The handy size, the extreme quietness of operation and the light sensor regulating the brightness of the screen are very pleasant for the patient as well as those in the vicinity.

Medela Chest Drainage University

Medela University is an online education platform that includes free-of-charge courses on chest drain management and in-service trainings, complemented by a variety of resources such as clinical competencies and the Thopaz+ Simulator.

REGISTER TODAY TO BUILD-UP AND REFRESH YOUR KNOWLEDGE.

Learn anytime, anywhere!

Find your topics of interest, add your selected courses to the shopping cart and you are ready for check out. All completed courses will be rewarded with a certificate.

<https://www.chestdrainage.medela.com>

1 Warnings and safety instructions

WARNINGS

Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTIONS

Indicates a potentially hazardous situation which, if not avoided, could result in minor or moderate injury.

Safety related tip

Indicating useful information about the safe use of the device.

Thopaz+ is approved exclusively for the use as described in these instructions for use. Medela can only guarantee the safe functioning of the system when Thopaz+ is used in combination with the original Thopaz+ accessories (canisters, tubings, carrying strap, mains adapter, docking station – see chapter 16 – Accessories/Appendix A).

Additional equipment connected to medical electrical equipment must comply with the respective IEC or ISO standards (eg. IEC 60950-1/ IEC 62368-1 for information and communication technology equipment). Furthermore all configurations shall comply with the requirements for medical electrical systems (see clause 16 of the IEC 60601-1). Anybody connecting additional equipment to medical electrical equipment configures a medical system and is therefore responsible that the system complies with the requirements for medical electrical systems. Attention is drawn to the fact that local laws take priority over the above mentioned requirements. If in doubt, consult your local Medela customer service representative or your technical service department.

Please read and observe these warning and safety instructions before operation. Please also familiarize yourself with associated information signals and troubleshooting instructions before operation (see chapter 12 – Troubleshooting and information signals). These instructions for use must be kept with the device for later reference.

Please note that these Instructions for use are a general guide for the use of the product. Medical matters must be addressed by a physician.

Medela considers itself only responsible for the effect on BASIC SAFETY, reliability and performance of Thopaz+ if it is used in accordance with the Instructions for use.

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

Any serious incident that has occurred in relation to the device should be reported to Medela AG and the relevant Competent Authority.

Subject to change.

WARNINGS

Warning: To reduce the risk of potential cross-contamination or exposure to biological hazards

- Wear gloves for all operations.
- After each use, the parts that have been in contact with aspirated secretions are to be cleaned and disinfected or disposed.
- The Thopaz⁺ system is not designed for retransfusion.

Warning: To reduce the risk of potential injury due to incorrect use

- Consult this manual prior to use.
- For use only by medically trained and qualified persons who have been adequately trained in the use of Thopaz⁺.
- The use of Thopaz⁺ for any other indication than intended is neither desired nor allowed.
- The pressure range to be set must be determined by a physician, based on the patient's age, weight, and clinical condition.
- For pediatric patients adapt pressure settings according to hospital guidelines.
- Note that smaller catheters may require higher pressure settings to ensure effective removal of fluids and bodily secretions.
- If a persistent air leak is shown, check if the system is assembled correctly before taking further corrective actions. Ensure the system is air tight by clamping the catheter and observing the air leak decrease to zero.
- Do not use Thopaz⁺ if the drainage therapy indicates a pressure greater than max. pressure range of -10 kPa.
- Do not use Thopaz⁺ if the drainage therapy indicates a flow rate greater than max. flow capacity of 5 l/min.
- Do not use Thopaz⁺ if the drainage therapy indicates no pressure should be applied to the patient.

Warning: To reduce the risk of potential injury during setup or operation

- "Pressure" generally implies "negative pressure".
- Thopaz⁺ was validated in combination with the accessories listed in Appendix A. For a correct and safe operation use Thopaz⁺ with these accessories only. Further information is supplied with the individual accessory.
- Do not connect bilateral thoracic drains to one Thopaz⁺ unit. In such cases, the use of two Thopaz⁺ units is recommended.
- To ensure Thopaz⁺ can reduce the pressure, it is essential to wait 30 seconds between taking a drainage sample and unclamping the patient catheter (drain).
- The canister is replaced on the basis of a visual check or according to the instructions on the Thopaz⁺ display (information signal).
- Thopaz⁺ system is regulated and controlled between the pump and the patient connector. Blockage and leakage at the catheter may not be detected by the Thopaz⁺ system, see chapter Troubleshooting and information signals. Regular patient monitoring is required by healthcare professionals.
- Leakage warning may not be triggered with set pressure below 0.4 kPa.

- Leakage warning is not triggered in case of disconnection of one catheter at the small double connector (079.0024).
- The catheter/connector interface is a location where clotting may occur. We recommend regular monitoring of this interface and appropriate removal procedure in case of occurrence.

Warning: To reduce the risk of potential injury due to interference with other devices

- Do not use Thopaz+ in MRT (Magnetic Resonance Tomography).
- HF (high-frequency) surgical equipment, radio networks or the like can influence the operation of the device and may not be operated in combination with the Thopaz+.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walki-talkie can affect the Thopaz+ and should be kept at least 30 cm away from the equipment.
- Data transfer via USB is not allowed during therapy.

Warning: To reduce the risk of potential electric shock or exposure to heat, fire or explosion

- Before cleaning the device, disconnect the device from any electrical plugs or connections.
- The device is not suitable for use while bathing, showering or in a hazardous explosive environment.
- Do not touch the patient and the docking contacts simultaneously.
- This product contains lithium-ion batteries which bear risk of fire, explosion and burns. Do not disassemble, crash, heat above 100 °C (212 °F), incinerate or dispose of in fire.
- Do not dry Thopaz+ with microwaves.

CAUTIONS

Caution: To reduce the risk of potential cross-contamination or exposure to biological hazards

- Visually inspect sterile packaging of the device for damage before opening.
Devices with a damaged packaging system must not be used
- Non sterile and reusable accessories must be cleaned and disinfected according to chapter 13 – General Reprocessing Guidelines.

Caution: To reduce the risk of potential injury due to incorrect use

- Incorrect use of Thopaz+ can cause pain and injury to the patient.
- The Thopaz+ 2l canisters are not intended for portable use (not by hand and not using the carrying strap).
- For safety reasons canister changes within the appropriate fluid level shall be recorded.
- Before diagnosing control detected fluid value for plausibility.
- For patients for whom a breakdown of the device can lead to a critical situation a replacement device must always be available.

Caution: To reduce the risk of potential injury during setup or operation

- No modification of this equipment is allowed.
- Before connecting Thopaz+ to the power supply, please verify that the supply voltage corresponds to that given on the device specification plate.
- For the initial startup do not switch Thopaz+ on before first charging the battery.

Safety related tips

- The product is delivered non-sterile. Before first use, as well as after each use, the product must be cleaned and disinfected according to a validated procedure.
- Never place Thopaz+ with an attached 2l canister on the docking station without attaching the optional adapter docking station 079.0038 for Thopaz+ 2l canister first.
- The acoustic warning “Leak in system – 301” is deactivated for approximately 6 minutes after Thopaz+ is switched on.
- Every time a new patient is connected, it is recommended to carry out the functional check.
- Before the first application in accordance with the intended use, an inspection according to IEC 62353 (see Appendix B) is recommended to acquire the reference values.
- Use sample port of tubing to drain a sample. Therefore refer to chapter 11 – Taking a drainage sample.
- Medela recommends to use the Thopaz tubing that best fits to the catheter.
- Make sure the sample port faces the patient.
- If the battery is low, use Thopaz+ with mains supply.
- The simethicone capsule requires 60 min of submersion for full effect.
- For most accurate fluid measurement place or keep Thopaz+ on the docking station (079.0037).

Safety instructions

- Thopaz+ is a medical device that requires specific safety measures in regard to EMC. It must be installed and put into operation in accordance with the attached EMC information in chapter 19 – Technical documentation.
- In each of the following cases, Thopaz+ must not be used and it must be repaired by Customer Services:
 - If the power cord or the plug are damaged
 - If the device is not functioning according to routine check
 - If the device is damaged
 - If the device shows clear safety defects.
- Thopaz+ has no user serviceable parts inside. For safety reasons, it is required that Thopaz+ is repaired throughout its service life strictly and exclusively by Medela authorised service centres.
- Keep the power supply cord away from hot surfaces.
- The mains plug and the on/off switch must not come into contact with moisture. Never pull the mains plug out of the fixed mains socket by pulling on the power supply cord.
- Thopaz+ must stand upright during use.
- Separation from the mains is only assured through the disconnection of the mains adapter and fixed socket connection.
- Never use the device at high room temperatures, while bathing or showering, if you are very tired or in an environment where there is a risk of explosion.
- Never place Thopaz+ in water or other liquids.
- When attaching a new canister, verify that the canister size displayed on screen is the same size as the canister attached.
- When using single use, sterile products, please note that they are not intended to be reprocessed. Reprocessing could cause loss of mechanical, chemical and/or biological characteristics. Reuse could cause cross contamination.
- Contact your local Medela customer service representative for assistance with product operations.
- The patient should be regularly monitored according to internal hospital guidelines.
- Use Thopaz tubing extension for adults only.

These instructions for use must be kept for later reference.

2 Power supply and battery operation

WARNINGS

This product contains lithium-ion batteries which bear risk of fire, explosion and burns. Do not disassemble, crash, heat above 100 °C (212 °F), incinerate or dispose of in fire.

Thopaz+ can be operated from the mains power supply [REF] 087.0059, from the docking station [REF] 079.0037 or with the integrated lithium-ion battery. Before using the docking station, read first the docking station instructions [REF] 200.1554. The battery is charged during mains operation. The battery duration is dependent on the run-time of Thopaz+. This is influenced by the extent of parenchymal leakage and the set pressure. Thopaz+ does not run continuously but only switches on when the actual and nominal values differ.

During continuous operation, Medela guarantees a minimum run-time of 4 hours in battery operation after the device has been fully charged. In practical operation the real battery run time results in >10 hours.

If Thopaz+ is not used very often, the battery must be charged approximately every 6 months to ensure optimum functioning.

Backup battery

If an internal defect occurs (broken cable, battery defective) Thopaz+ turns off and the acoustic warning is sounded for at least 3 minutes (powered by a backup battery). Under these conditions Thopaz+ functions as a one way valve. Replace Thopaz+ immediately.

Charge battery

Plug the Thopaz+ power supply into a wall outlet and open the small protection cover to access the DC socket. Plug the coaxial power supply into the DC socket. Moving bars on the battery display will indicate charging is occurring [].

Battery low


Approximately 30 minutes before the battery is fully empty an acoustic signal sounds and the battery symbol [] starts to blink. The acoustic signal can be muted but the battery symbol continues to blink until the battery is empty. The set pressure is maintained but the battery should be recharged as soon as possible.

Press simultaneously to acknowledge an acoustic information signal and follow instructions on the display (see chapter 12 – Troubleshooting and information signals).

Battery empty

The battery symbol will blink with an acoustical sound for 10 minutes prior to the battery fully discharging, unless the pump is switched off earlier. The acoustic sound cannot be muted during the 10 minutes. If Thopaz+ is not connected to a power source the pump will switch off automatically after 10 minutes. The negative pressure is not maintained after the pump turns off.

1. Battery empty

2. Battery almost empty

3. Battery is being charged (moving bars)

4. Battery is fully charged and connected to the mains

* Symbol blinks

3 Description

Introduction

Thopaz+ is a high-quality digital chest drainage system. The compact system provides regulated negative pressure close to the patient's chest. It monitors the air leak and the fluid collection. Thopaz+ only applies vacuum required to maintain the negative pressure prescribed and set by the managing physician. The digital colour display provides objective data in real time as well as in historical graphs, which allows easy tracking of the therapy progress. It combines easy handling and cleaning with safety features to ensure optimal operation.

Intended use/purpose

Thopaz+ is intended to be used for aspiration and removal of surgical fluids, tissue, gases, bodily fluids or infectious materials.

Indications for use

Thopaz+ is indicated for all situations where chest drains are applied – especially for thoracic drainage in the pleural and mediastinal cavity in situations such as pneumothorax, after cardiac or thoracic surgery (post-operative), thorax injury, pleural effusion, pleural empyema or other related conditions in which regulated pressure is indicated. Thopaz+ is intended for use on patients in appropriate care settings.

Contraindications

There are no known contraindications for the Thopaz+ digital chest drainage system.

Intended user

Thopaz+ should only be operated by healthcare professionals familiar with cardiothoracic drainage. These persons must not be hard of hearing or deaf and must have adequate visual faculty.

Intended patient population

Thopaz+ is intended to be used on patients only exhibiting conditions as described in the indications for use.

Important note

Compliance with proper surgical procedures and techniques is the responsibility of the physician. Each physician must evaluate the appropriateness of the treatment based on his own knowledge and experience.

Undesirable side effects



There are no known undesirable side effects associated with the Thopaz+ digital chest drainage system. However, the following complications may occur during cardiothoracic drainage:

- Pleural effusion (may lead to e.g. chest infection)
- Extra-pulmonary accumulation of air (may lead to e.g. pneumothorax/tension pneumothorax, subcutaneous emphysema)
- Pericardial accumulation of fluid (may lead to e.g. cardiac tamponade)
- Pleurocutaneous fistula
- Re-expansion pulmonary oedema
- Chest pain/discomfort
- Coughing


4 Overview

⚠ WARNINGS
Do not touch the patient and the docking contacts simultaneously.

Main elements of Thopaz+

Display

How to activate a button


Press the buttons briefly to activate the desired action.

Exception: for the function switch to “Standby”, press the button for longer than 3 seconds.

Tubings

STERILE EO

Normal/small/large/single/double connector
 Material: PVC (medical grade) sterile, twofold packaging
 Length: 1.5 m

079.0025/26: Patient tube inner diameter: 5.5 mm
 Connector outside diameter 1: 9.5 mm
 Connector outside diameter 2: 16.0 mm

079.0021/22: Patient tube inner diameter: 5.5 mm
 Connector outside diameter 1: 6.4 mm
 Connector outside diameter 2: 16.0 mm

079.0023/24: Patient tube inner diameter: 5.5 mm
 Connector outside diameter 1: 4.0 mm
 Connector outside diameter 2: 11.0 mm

Thopaz tubing single,
large connector
REF 079.0025

Thopaz tubing single
REF 079.0021

Thopaz tubing single,
small connector
REF 079.0023

Thopaz tubing double,
large connector
REF 079.0026

Thopaz tubing double
REF 079.0022


Thopaz tubing double,
small connector
REF 079.0024*

*Leakage warning is not triggered in case of disconnection of one catheter at the small double connector (079.0024).

Safety related tip

Medela recommends to use the Thopaz tubing that best fits to the catheter.
Make sure the sample port faces the patient.

Sample port

Place tubing with sample port
facing the patient.

Use the Thopaz tubing that best
fits to the catheter.

Canisters


STERILE EO

With/without solidifier 0.3l/0.8l/2l

Material: Polypropylene, sterile

Solidifier: 0.3l = 9g/0.8l = 26g/2l = 2 x 30g

Simethicone capsule: Supports reliability of digital fluid measurement (0.3l = 720mg/0.8l = 1440mg/2l = 1440mg Simethicone)

Thopaz canister 0.8l, antifoaming

REF 079.0016

Thopaz canister 0.8l with solidifier, antifoaming

REF 079.0017

Thopaz canister 0.3l, antifoaming

REF 079.0011

Thopaz canister 0.3l with solidifier, antifoaming

REF 079.0012

Thopaz canister 2l, antifoaming

REF 079.0018

Thopaz canister 2l with solidifier, antifoaming

REF 079.0019

Safety related tip

Use sample port of tubing to drain a sample.

Therefore refer to chapter 11 – Taking a drainage sample.

Safety related tip

The simethicone capsule requires 60 min of submersion for full effect.

5 Additional information

Definition of vacuum

By the application of medical aspiration devices, vacuum is normally given as the difference (in absolute figures) between absolute pressure and atmospheric pressure or as negative values in kPa. In this document, the indication of -10 kPa for example always refers to a pressure range in kPa below atmospheric ambient pressure (according to EN ISO 10079:1999).

Regulated pressure

The pressure is regularly checked and adjusted by Thopaz⁺. The set pressure corresponds to the pressure of the entire Thopaz⁺ system.

Safety pressure limit

The adjustable pressure range for Thopaz⁺ is $-0,1$ to -10 kPa (-100 mbar, -75 mmHg, -102 cmH₂O). If pressure values greater than -7 kPa are set, the following warning appears in the display: "too high pressure levels can cause pain and/or serious injuries to the patient". This warning must be acknowledged with "OK" before the pressure can be increased further.

Tube flushing

Flushing of the tubing occurs every 5 minutes or when Thopaz⁺ detects a siphon and prevents clogging of the patient tubing.


Safety chamber in canister

For appropriate, safe operation Thopaz⁺ must be in an upright position. If Thopaz⁺ tips over, the special construction of the safety chamber in the upper region of the canister protects the hydrophilic filter from premature clogging. The pressure is maintained. Important: immediately set Thopaz⁺ upright again. This function is disabled if the maximum filling capacity of the canister is exceeded.

Fluid measurement

Thopaz+ calculates the fluid amount from the fluid level in the canister and in regard to the attached canister size. The fluid level in the canister gets detected by capacitive coupling. For a valid fluid measurement Thopaz+ needs to stay within a tolerance of 10° vertical (upright). Fluid measurement starts a latest with a minimal fluid amount of:

Thopaz canister without solidifier	Thopaz canister with solidifier
30 ml in a 0.3l canister	40 ml in a 0.3l canister
50 ml in a 0.8l canister	70 ml in a 0.8l canister
200 ml in a 2.0l canister	

If these conditions are not fulfilled the fluid values appear in grey with a crossed out fluid symbol []. Additionally a Thopaz+ out of axis symbol [] or a Thopaz+ not enough fluid symbol [] appears on the action bar. No additional fluid can be measured in this status. The grey fluid values show the last measured fluid values.

This type of fluid measurement allows a measurement accuracy of ± 15 ml in a 0.3l canister, ± 40 ml in a 0.8l canister and ± 100 ml in a 2.0l canister (corresponding to $\pm 5\%$ of the canister volume).

Thopaz+ displays the fluid values in a 0.3l canister in 1 ml steps, in a 0.8l canister in 5 ml steps and in a 2.0l canister in 10 ml steps. After a total drainage amount of 5000 ml and more the display will switch to liter (l). Thopaz+ will display fluid values in 0.01 l steps (equals 10 ml steps) thereupon.

Limitations

Thopaz+ calculates the fluid measurement according the fluid level in the respective canister. Movement of Thopaz+ can cause fluids to temporary cover the canister walls. A potential fluid film covering the canister wall can result in temporary higher fluid volume read out displayed on the screen. Also foam formation in the canister can result in higher fluid volume read out displayed on the screen. All Thopaz canisters (079.0011; 079.0012; 079.0016; 079.0017; 079.0018; 079.0019) are equipped with one or two simethicone capsules containing 720 mg simethicone each (see chapter 4 – Overview). Simethicone prevents foam formation and the coverage of fluid films on the canister wall.

Safety related tip

The simethicone capsule requires 60 min of submersion for full effect.

Safety related tip

For most accurate fluid measurement place or keep Thopaz+ on the docking station (079.0037).

Overflow protection/bacteria filter

A hydrophilic overflow protection/bacteria filter in the canister as well as in the tube sets (in the connector to Thopaz+) protects Thopaz+ from the entry of fluids and protects the environment from contamination.

Light sensor

Thopaz+ is equipped with a light sensor. It adjusts the display brightness to the surrounding light conditions. Thus, the patient does not get disturbed while sleeping by a bright shining display.

6 Installation

6.1 Check initial delivery

Check the delivery package of Thopaz+ for completeness and general condition.

Thopaz+

REF 079.1000/101046671* with Instructions for use 200.6841
 079.1002 with Instructions for use 200.6842
 079.1003 with Instructions for use 200.6843/101034256
 *not CE marked

Thopaz+/Liberty wall mains adapter/charger

REF 087.0059

Instructions for use

REF 200.6841 (for 079.1000/101046671)
 200.6842 (for 079.1002)
 200.6843/101034256 (for 079.1003)

6.2 Initial startup

WARNINGS

For use only by medically trained and qualified persons who have been adequately trained in the use of Thopaz+.

Safety related tip

Before the first application in accordance with the intended use, an inspection according to IEC 62353 (see Appendix B) is recommended to acquire the reference values.

Important: Do not connect the patient to the system yet!

1.

- 1.1 Select plug.
- 1.2 Attach plug to the wall mains adapter
REF 087.0059.

2.

- 2.1 Connect Thopaz+ to mains by mains adapter
REF 087.0059 or
by docking station
REF 079.0037.

3.

- 3.1 Charge for approximately 3 hours or until the [] symbol appears in the display.
- 3.2 Press [] to switch Thopaz+ on. A selftest starts.

4.


- 4.1 Select language.
- 4.2 Confirm with "OK".
- 4.3 Follow instructions on display.

5. 5.1 Press [] to switch Thopaz+ off.

6. 6.1 Disconnect Thopaz+ from the fixed mains socket by pulling on the plug housing. Do not pull on the cable or the bend protection.

7. Thopaz+ is ready for preparation for use.

7 Preparation for use/operating instructions

WARNINGS

- For use only by medically trained and qualified persons who have been adequately trained in the use of Thopaz+. Wear gloves for all operations. “Pressure” generally implies “negative pressure”.
- Adapt pressure settings for pediatric patients.
- Note that smaller catheters may require higher pressure settings to ensure effective removal of fluids and bodily secretions.

CAUTIONS


- Visually inspect sterile packaging of the device for damage before opening. Devices with a damaged packaging system must not be used.
- Non sterile and reusable accessories must be cleaned and disinfected according to chapter 13 – General Reprocessing Guidelines.
- For the initial startup do not switch Thopaz+ on before first charging the battery.

7.1 Checks before use

- Check the Thopaz+ system before use for damage of the power cord or plug, obvious device damage or safety defects and proper functioning of the device.
- Check that the seal is correctly placed in the vacuum port of Thopaz+. If you do not know where the vacuum port is placed please refer to chapter 4 – Overview.
- Make sure that the rechargeable battery is recharged in case Thopaz+ is operated in battery mode.
- Check all accessories prior to use:
 - Canister for cracks, brittle and flawed spots.
 - Tubing for cracks, brittle areas and that connectors are firmly attached.Replace if necessary.

7.2 Connect tubing

1.

- 1.1 Open external packaging.
- 1.2 Keep patient connector in internal bag for hygienic reasons.

2.

- 2.1 Insert the smaller of the two connectors (connection to pump) horizontally in the direction of the arrow.

7.3 Snap in canister


- 1.1 Unpack the canister.
- 1.2 Position the openings at the top and attach the bottom part of the canister to Thopaz+.

- 2.1 Push the canister into Thopaz+ until you hear a click.

7.3.1 Using 2l canister

- 1.1 Unpack the canister.
- 1.2 Position the openings at the top and attach the bottom part of the canister to Thopaz+.
- 1.3 Push the canister into Thopaz+ until you hear a click.

Safety related tip

Never place Thopaz+ with an attached 2l canister on the docking station without attaching the optional adapter docking station **REF** 079.0038 for Thopaz+ canister 2l first.

CAUTIONS


The Thopaz+ 2l canisters are not intended for portable use (not by hand and not using the carrying strap).

7.4 Switch Thopaz+ on

Important: Do switch Thopaz+ on before connecting the system to the patient.

1.  1.1 Press [] to switch Thopaz+ on. A selftest starts.
1.2 Check acoustic signal (beep).


If the selftest is not successful, follow the troubleshooting instructions on the display and refer to chapter 12 – Troubleshooting and information signals.

2.  2.1 Confirm whether a new patient is connected or not. The patient number is important for therapy time, history data and data transfer to the PC.


yes Thopaz+ will issue a new patient number. The therapy time and history data will start from 0 (recommended for a new patient).

Note: Once the patient number reaches 100 it automatically defaults back to 1.

no The patient number remains unchanged (recommended for the continued treatment of the same patient). The data history and the therapy time are retained.

3.  3.1 Check canister matching: detected canister size vs. attached canister. No match: dispose canister.

7.5 Carry out functional check

1.1 Seal off the conical connector of the patient tubing with your thumb (through the packaging). The tubing must not be kinked otherwise the measuring tube will be closed.

2.1 Press "OK".

Test passed:

Thopaz+ is ready for use.

Test failed:

Leak in the system:

- Check tubing connection.
- Check seal is correctly placed.
- Check pressure relief valve for tight connection.
- Replace canister and/or tubing if necessary.

3. Thopaz+ is in Standby and is ready for use with the factory settings. (To change factory settings please refer to chapter 10 – Changing settings)
4. Connect the Thopaz+ system to the patient catheter (drain) according to hospital guidelines. To start the therapy press "On".

Safety related tip

The acoustic warning "Leak in system – 301" is deactivated for approximately 6 minutes after Thopaz+ is switched on.

Safety related tip

Every time a new patient is connected, it is recommended to carry out the functional check.

7.6 Check therapy progress

7.6.1 Main display

The current size of the air leak (ml/min) is shown on the left side. The total amount of collected fluid (ml) during the therapy is shown on the right.

According to your preference the air leak history of 12h/100 ml/min and two customized fluid parameters can be displayed in addition. (To change factory settings please refer to chapter 10 – Changing settings).

- Air leaks between 0 and 1000 ml/min are displayed in 10 ml steps.
- Air leaks over 1000 ml/min are displayed in 100 ml steps.
- The air leak is measured with a tolerance of $\pm 20\%$.
- Rounding function

The digital display of air leak values is rounded to zero ml/min (0 ml/min) if the air leak measurement indicates an air leak < 5.0 ml/min. If the air leak measurement indicates an air leak ≥ 5.0 ml/min and < 15 ml/min the digital display of air leak values is rounded to ten ml/min (10 ml/min)

- The fluid measurements are displayed as follow: In a 0.3 l canister in 1 ml steps, in a 0.8 l canister in 5 ml steps and in a 2.0 l in 10 ml steps. After a total drainage amount of 5000 ml and more the display will switch to liter (l). Thopaz+ will display fluid values in 0.01 l steps (equals 10 ml steps) thereupon.

7.6.2 Reset fluid display

For an individual fluid management an additional fluid parameter can be set to zero at any given time.

1. 
 - 1.1 Press [] for longer than 3 seconds.

2. 
 - 2.1 The timer on the left side is giving you the time when the parameter was set to zero (in this case 12 minutes ago). The fluid amount on the right side is giving you the drained fluid volume (in this case 50 ml got drained in the last 12 minutes). To reset the timer and fluid again press [] for longer than 3 seconds.

7.6.3 History


- Thopaz+ provides air leak history graphs and fluid history graphs.
- For scrolling through the graphs press "Next".
- To get back to your main display press "Home".

Air leak history

Press [] to get into the air leak history.

- Three different air leak history graphs (72 h/Autoscale, 24 h/1000 ml/min, 24 h/100 ml/min) plus the catheter check (see chapter 7.7 – Catheter check) are available.
- The air leak autoscale history graph adjusts the air leak scale to the largest measured air leak from the last 72 h to 1000/2000/3000/4000/5000 ml/min.

Pressure scale

Air leak scale (in ml/min)

Air leak

Current time

Current page



Set and measured pressure
(depending on your settings)

Fluid history

Press [] to get into the fluid history.

- Three different fluid history graphs are available. 72 h/Autoscale, 24 h/100 ml and 6 h/Autoscale.
- The fluid autoscale history graph shows the total collected fluid amount over the last 72 h.


7.7 Catheter check

Important: the catheter check works only with a patient catheter (drain) in the interpleural space.

Press [X] to get into the air leak history. Scroll to page 4/4 by pressing "Next" [Next].

- The catheter check is only active if the air leak is 0 ml/min.
- Follow instructions on screen.
- A pressure difference from minimal 3 cmH₂O between inspiration [▲] and expiration [▼] is displayed at an air leak of 0 ml/min. This pressure difference is an indication for the respiration cycle of the patient and confirms that the patient catheter (drain) is not clogged.

▲ Patient catheter (drain) open

▼ Patient catheter (drain) clogged?


This number shows the measured maximal and minimal pressure of each breathing cycle. This pressure is indicated in the set pressure unit.

7.8 Switching Thopaz⁺ off

1. 
 - 1.1 Clamp patient tubing with the tube clamp.
 - 1.2 Clamp patient catheter (drain).
 - 1.3 Press "Standby" for longer than 3 seconds to switch off the pressure.
2. 
 - 2.1 Press [] to switch off Thopaz⁺.
3. Release, remove and seal the canister by the canister seal. Dispose of canister and patient tubing in accordance with internal hospital guidelines.
4. Clean Thopaz⁺ according to chapter 13 – General Reprocessing Guidelines.

7.9 Transferring data to PC with ThopEasy+

With ThopEasy+ all stored data can be transferred to a PC as documentation, and for completion of the patient records. The data can be expanded with patient information, saved and printed. Operating instructions and the software can be found on the Medela website.

For the connection between Thopaz+ and a PC, use the USB cable [REF 079.0034](#).

ThopEasy+ is compatible with the following Microsoft Operating Systems:

- Windows 7 32 Bit
- Windows 7 64 Bit
- Windows XP 32 Bit (SP3 or higher)
- Windows XP 64 Bit (SP3 or higher)
- Windows 8 32 Bit
- Windows 8 64 Bit
- Windows 8.1 32 Bit
- Windows 8.1 64 Bit
- Windows 10 64 Bit
- Windows 11 64 Bit

WARNINGS

Data transfer via USB is not allowed during therapy.

8 Replace canister

WARNINGS

The canister is replaced on the basis of a visual check or according to the instructions on the Thopaz+ display (warning signal). Wear gloves for all operations.


CAUTIONS

Canister changes shall be recorded for safety reasons.

1. Prepare sterile canister (with/without solidifier).
2. Clamp off patient tubing with the tube clamp.

- 3.1 Switch Thopaz+ to Standby by pressing "Standby" for longer than 3 seconds.

- 4.1 Release the canister by pressing the release button.
- 4.2 Remove canister.

- 5.1 Unpack the canister.
- 5.2 Position the openings at the top and attach the bottom part of the canister to Thopaz+.

- 6.1 Push the canister into Thopaz+ until you hear a click.

Safety related tip

Using a 2l canister refer to chapter 7.3.1 – Using 2l canister.


7.1 Check canister matching:
Detected canister size vs.
attached canister.
No match: Dispose
canister.

7.2 Press "On". Pressure is
built up and unclamp the
tube clamp.

8.1 Check credibility of air
leak value.


9.1 Seal the used canister with
the canister seal.

10. Dispose of the used canister in accordance with internal
hospital guidelines.

Canister with solidifier: activation of the solidification process

1.1 Check if canister is sealed
with the canister seal.
1.2 Press – to open the
solidifier chamber.

2.1 Shake – to activate the
solidification process.
2.2 Dispose of the used
canister in accordance
with internal hospital
guidelines.

9 Changing pressure during operation

WARNINGS

The negative pressure may only be changed by the physician or on medical orders.

Change pressure

Important: Thopaz⁺ is running.

Press simultaneously

- 1.1 Press [▶] and [◀] simultaneously.
- 1.2 Change pressure by pressing [↕] or [↔] and confirm with "OK" [◀].

9.1 Physiological pressure

For patients who are to be treated by gravity drainage (= water seal), the physiological pressure can be activated.

- 1.1 Press "Physio" [▶].
- 1.2 Confirm with "OK" [◀].

This mode corresponds to a pressure of -0.8 kPa / -6 mmHg / $-8 \text{ cmH}_2\text{O}$ / -8 mbar .

10 Changing settings

10.1 Factory settings

Thopaz+ is delivered with the following factory settings. These settings can be changed and saved as the new standard.

- Pressure:** -2.0 kPa
- Choose to max. -10 kPa (-102 cmH₂O, -100 mbar, -75 mmHg).
- Pressure unit:** kPa
- Choose between kPa, cmH₂O, mmHg, mbar.
- Measured pressure:** Hide
- Choose between hide and show (only relevant in the air leak history).
- This feature is added for advanced Thopaz+ users. A higher measured pressure than set pressure can be an indication that the patient is able to regulate the negative pressure in the inter-pleural space.
- Fluid time 1:** OFF
- Choose between 24 h/12 h/8 h/4 h/OFF
- Shows the total amount of collected fluid over the last X hours.
- Fluid time 2:** OFF
- Choose between 2 h/1 h/0.5 h (30 min)/0.25 h (15 min)/OFF
- Shows the total amount of collected fluid over the last X hours/minutes.
- Fluid alarm:** OFF
- Choose between 500 ml/5 min, 15 min or 60 min; 400 ml/5 min, 15 min or 60 min; 300 ml/5 min, 15 min or 60 min; 200 ml/5 min, 15 min or 60 min; 100 ml/5 min, 15 min or 60 min; 50 ml/5 min, 15 min or 60 min/OFF
- The fluid alarm signal 407 appears if Thopaz+ detected a higher fluid collection in the last 5 min, 15 min or 60 min than the desired parameter.
- Air leak history:** Hide
- Choose between hide and show.
- This history can be displayed on the main display (see chapter 7.6.1 – Main display).
- Language:** English
- Choose your language.
- Save as standard:** Choose if you want to save the new parameter as a standard (guidance see chapter 10.3 – Save a changed parameter as the new standard).

10.2 Change settings

WARNINGS



The settings may only be changed by the physician or on medical orders.

CAUTIONS

Adapt pressure settings for pediatric patients.

Important: Thopaz+ is on Standby.

1.  1.1 Press "Menu" [].
2.  2.1 Choose the desired parameter by pressing [].
2.2 Confirm with "OK" [] to change the desired parameter.
2.3 Change the desired parameter by pressing [] or [] and confirm with "OK" [].
2.4 Press "Home" [] to leave the settings.
2.5 The settings are changed for the current therapy.

10.3 Save a changed parameter as the new standard


1. Change all desired parameters as mentioned above.
2. Choose "Save as standard" and press "OK" []. "Saved" will appear for approx. 3 seconds.
3. The parameters are saved as the new standard. Whenever you switch Thopaz+ on and choose a new patient, Thopaz+ will be set with these settings.
4. Press "Home" [] to leave the settings.

11 Taking a drainage sample

WARNINGS

To ensure Thopaz+ can reduce the pressure, it is essential to wait 30 seconds between taking the sample from the sample port and unclamping the patient catheter (drain).


1. Before taking a drainage sample make sure that there is fluid in the patient tube.

2. Clamp the catheter (drain).

3.1 Switch Thopaz+ to Standby by pressing "Standby" for longer than 3 seconds.

4. Disinfect the sample port before taking the sample with CaviWipes or Mikrozid AF Wipes.

5.1 Remove with syringe* air from the patient tubing. Repeat until fluid has gathered at the sample port.

* Syringe 17G (1,4mm) or thinner.

6. With a syringe* take a sample from the patient tube.

7.1 Switch Thopaz+ on by pressing "On" – pressure is built up.

8. **Wait 30seconds!**

To ensure Thopaz+ can reduce the pressure, it is essential to wait 30seconds between taking the sample and unclamping the patient catheter (drain).

9. Unclamp the patient catheter (drain).


12 Troubleshooting and information signals

Description

Thopaz+ distinguishes between cautions (yellow), warnings (red) and internal errors (red). If Thopaz+ detects any of these situations, an acoustic signal sounds and a description of the problem appears on the display. By pressing the two selection buttons simultaneously, the acoustic information signals are suppressed for 60 seconds. Exceptions: For the caution 306 "canister full" and the warnings 302 "tubing clogged", 313 "filter clogged" and 315 "Thopaz+ overheated" the acoustic information signals are suppressed for 5 minutes. Information signals do not replace the need of regular monitoring of the patient. Leakage warnings may not be triggered and therefore the device does not alert in case of disconnection of the Thopaz tubing in the following situations

- set pressures below 0.4 kPa
- disconnection with the small double connector (079.0024)

Example:

Visible on display					
Error number	Problem description	Troubleshooting instructions	Remarks / Potential source of error	Pressure	
Information signal: Caution (yellow)	306	Canister full	1. Clamp tubing 2. Switch to Standby 3. Replace canister 4. Continue with "On" 5. Unclamp tubing		Yes
	401	Battery low	Connect Thopaz+ to mains	Remaining time of battery approx. 30min.	Yes
	402	USB connection not permitted	Unplug USB cable	USB connection not permitted during operation or in the vicinity of the patient.	Yes
	405	Standby	Switch Thopaz+ on or off	In Standby for 5minutes.	No
	406	Temperature too high	Prepare backup device	Do not expose Thopaz+ to heat sources.	Yes
	409	Fluid level sensor out of order	– Fluid measurement unavailable – Confirm and continue therapy		Yes
	410	Canister detection	1. Clamp tubing 2. Switch to Standby 3. Replace canister 4. Continue with "On" 5. Unclamp tubing		Yes
Information signal: Warning (red)	301	Leak in system	Check system for leaks	– Disconnection. – Missing seal at vacuum port. – Thopaz+ detects leaking in the Thopaz+ system including Thopaz tubing, NOT of disconnection of the patient catheter from the patient. – Leakage warning may not be triggered with set pressure below 0.4kPa. – Disconnection of one catheter at the small double connector will not be detected by Thopaz+	No
	302	Tubing clogged	1. Clamp catheter 2. Clamp tubing 3. Switch to Standby 4. Replace tubing 5. Continue with "On" 6. Unclamp catheter	– Tubing is kinked or clogged for approx. 10min. – Thopaz+ detects clotting of Thopaz+ tubing, NOT of the patient catheter.	No
	305	Battery empty	Connect Thopaz+ to mains		No

Visible on display					
	Error number	Problem description	Troubleshooting instructions	Remarks / Potential source of error	Pressure
Information signal: Warning (red)	311	Selftest failed	Snap the canister out and in again	Do not switch Thopaz+ on if the patient is already connected.	No
	313	Filter clogged	1. Clamp tubing 2. Replace canister 3. Continue with "On" 4. Unclamp tubing	Filter in canister clogged.	No
	315	Thopaz+ overheated	Replace Thopaz+		No
	407	Fluid alarm	Patient is losing a large amount of fluid		Yes
Internal error		Internal error XXX	1. Switch Thopaz+ off and on again 2. Inform Medela Customer Service	If error message repeats, inform Medela Customer Service.	

Safety related tip

The acoustic warning "Leak in system – 301" is deactivated for approximately 6 minutes after Thopaz+ is switched on.

Safety related tip

Caution "canister full"/Warning "filter clogged".

The caution "canister full" is triggered when the canister filling level approximately reaches the maximal scale line indicated on the respective canister label. Change the canister according the instructions of use, chapter "8 – Replace canister". If secretions line the walls of the canister, the caution "canister full" may be activated prematurely.

After muting the caution with the selection buttons, a symbol [] appears during operation. After 5 minutes the warning reappears acoustically.

An overfilling of the respective canister volume will trigger the warning "filter clogged". The warning is triggered due to the clogging of the hydrophilic filter upon contact with liquids or secretions. Please refer to chapter "4 – Overview" for location of the hydrophilic, overflow protection/bacteria filter in the canisters. Subsequent to a clogged filter and the warning "filter clogged" the vacuum applied to the patient is interrupted and cannot be regulated or maintained.

An additional accumulation of liquids, secretions and air is released through the pressure relieve valve at the positive pressure of 0.2–0.5 kPa. This pressure release function is not a substitute for a free overflow function. Therefore the canister must be replaced once full. Please refer to chapter "4 – Overview" for location of the pressure relieve valve.

Change the canister according to the instructions for use, chapter "8 – Replace canister".



13 General reprocessing guidelines

Thopaz ⁺ , docking station, standard holder, and holder – standard reprocessing				
Thopaz ⁺ in case of spill on seal				
In case of spill on release button, bed holder, docking station's push button standard holder, and holder with rail				
Additional information	x	x	x	<ul style="list-style-type: none"> – Per ISO 17664-2, these instructions have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed by using equipment, materials and personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process. – For specification of water qualities see AAMI TIR34.
	x	x	x	<ul style="list-style-type: none"> – The washer-disinfector shall be qualified according to ISO 15883 series. Cleaning and disinfection was validated in a washer-disinfector of the type: Uniclean PL II 1-2 EL (MMM, Germany). – All disassembled parts must be safely fixed in the carriers/fixation points. – Do not overload the washer-disinfector. Arrange the disassembled parts in such a way that no areas are left unwashed and inner and outer surfaces are reached by the cleaning liquids.
Limitations	x			<ul style="list-style-type: none"> – Discard or service the device (or component as applicable) if it shows visible signs of wear or damage.
General safety instructions	x	x	x	<ul style="list-style-type: none"> – Always wear personal protective equipment (PPE): disposable gloves and other PPE as per local guidelines and regulations. – Perform point-of-use treatment with utility water ($\leq 40^{\circ}\text{C}$, $\leq 104^{\circ}\text{F}$) directly after use of the device (before soil can dry onto the device). Violation of this may result in the fixation of residue and thus inhibit disinfection. – If the product is used on a patient who suffers from a disease, and whose pathogens cannot be eliminated with the procedure outlined below, the product must be disposed of. – Consult the cleaning and disinfection agent manufacturer's instruction for use regarding, including but not limited to exposure times and safety measures.
	x	x		<ul style="list-style-type: none"> – Never immerse the device in or rinse with water or other liquids. Do not spray cleaning agent and disinfectant directly on the device.

Thopaz⁺, docking station, standard holder, and holder – standard reprocessing

Thopaz⁺ in case of spill on seal

In case of spill on release button, bed holder, docking station's push button standard holder, and holder with rail

Preparation before cleaning	x	x	x	– Disconnect/remove all accessories.
	x	x		– Wipe external surfaces of the device to remove all gross soil with a soft, lint-free wipe moistened with utility water.
			x	<ul style="list-style-type: none"> – Disassemble all soiled components into their maximum number of individual parts. – Remove the release button by unlocking it with the help of two levers, that have no sharp edges and a shaft diameter of ≤ 3.5 mm (see picture). – Loosen the screws of the bed holder with a Philips screwdriver and remove them (see picture). <div style="display: flex; justify-content: space-around;">   </div> <ul style="list-style-type: none"> – Carefully open the Torx Screw on the docking station, compress the spring by pressing the button. After removing the screw, slowly release the push button. Next, remove the push button and the spring. Then remove the lower claw by closing the clamp and then pulling.
		x	<ul style="list-style-type: none"> – If necessary, and for the removal of gross soil, place the disassembled components in utility water for 10 minutes and wipe off visible staining with a soft, lint-free wipe soaked in utility water. – If residual soil has dried onto the disassembled components, the soil must be rehydrated before the enzymes can be effective. 	
Manual cleaning			x	<ul style="list-style-type: none"> – Clean thoroughly under the gasket with the help of a sterile spatula without sharp edges and according to the following process: – Fold a CaviWipes™ or Incidin OxyWipe S™ wipe once in the middle so that the tip of the spatula is enveloped by the wipe. – Wrap the wipe once around the spatula. – Carefully insert the wipe-wrapped spatula under the contaminated seal on the far left or far right side of the Thopaz⁺. Then slide it to the center of the device and push it outward in an sweeping motion. – Replace the cleaning and disinfection wipe with a new one. – Carefully insert the wrapped spatula under the contaminated seal on the other side of Thopaz⁺. Then push towards the center of the device, past the spot where you swept out before. Then slide outwards in a sweeping motion. – Replace the cleaning and disinfection wipe with a new one and repeat the aforementioned cleaning steps until the cleaning and disinfection wipe is no longer soiled.
		x	x	<ul style="list-style-type: none"> – Wipe all external surfaces of the device with CaviWipes™ or Incidin OxyWipe S™ wipe. – Wipe away from difficult-to-clean areas (e.g., where components that cannot be disassembled meet and crevices). – Use a new cleaning and disinfectant wipe when the wipe is contaminated. – Clean until all visible soil is removed. Pay special attention to the difficult-to-clean areas.

Thopaz+, docking station, standard holder, and holder – standard reprocessing

Thopaz+ in case of spill on seal


In case of spill on release button, bed holder, docking station's push button standard holder, and holder with rail

Manual disinfection	x	<ul style="list-style-type: none"> - Disinfect thoroughly under the gasket with the help of a sterile spatula without sharp edges and according to the following process: - Fold a new CaviWipes™ or Incidin OxyWipe S™ wipe once in the middle so that the tip of the spatula is enveloped by the wipe. - Wrap the wipe once around the spatula. - Carefully insert the wrapped spatula under the contaminated seal on the far left or far right side of the Thopaz+. Then slide it to the center of the device and push it outward in an arcuate motion. - Replace the cleaning and disinfection wipe with a new one. - Carefully insert the wrapped spatula under the contaminated seal on the other side of Thopaz+. Then push towards the middle of the device, past the spot where you swept out before. Then slide outwards in a sweeping motion.
	x x	<ul style="list-style-type: none"> - Take a new CaviWipes™ or Incidin OxyWipe S™ wipe and wipe all external surfaces of the equipment. - Pay special attention to the difficult-to-clean areas of the device. - To aid exposure of difficult-to-clean areas, a new CaviWipes or Incidin Oxy Wipe S wipe may be wrapped around a spatula or a similar utensil. - Make sure all surfaces of the device remain visibly moist at room temperature for 3 minutes. If the surface is getting too dry, moisten the surface using a new wipe. - Take a new CaviWipes™ or Incidin OxyWipe S™ wipe and wipe again all external surfaces of the equipment. - After the prescribed exposure, remove any residuals using a soft, lint-free wipe moistened with critical water per AAMI TIR 34.
Automated cleaning	x	<ul style="list-style-type: none"> - If applicable, place disassembled components in a small sieve tray, which is positioned on the load carrier. - Place the disassembled parts in a flat tray so that they cannot flip over. Orient with recesses facing down that liquids are drained by gravity. - Do not use any drying aids (rinsing agents). These could remain on the surface with a detrimental effect to the device and its biocompatibility. <p>The cleaning program in the washer-disinfector consists of the following:</p> <ul style="list-style-type: none"> - 1 minute pre-cleaning with utility water - 5 minutes cleaning at 55 °C with 0.5 % solution of neodisher® MediClean forte (alkaline enzymatic cleaner) in utility water - 1 minute rinsing with critical cold water per AAMI TIR 34
Automated disinfection	x	<ul style="list-style-type: none"> - Thermal disinfection with critical cold water per AAMI TIR 34 (without a chemical disinfectant) at 90 °C for 1 minute (A₀=600) or adapt A₀ values per local guidelines and regulations.

Thopaz⁺, docking station, standard holder, and holder – standard reprocessing

Thopaz⁺ in case of spill on seal

In case of spill on release button, bed holder, docking station's push button standard holder, and holder with rail

Drying		x	– Dry disassembled components in washer-disinfector at 110 °C for at least 45 minutes.	
	x	x	x	– If drying in the washer-disinfector is not possible or in case of residual moisture, wipe external surfaces dry using a dry, soft, lint-free wipe, or carefully dry with medical grade compressed air. – Pay special attention to the dryness of hard to reach areas.
Inspection	x	x	x	– Visually inspect the device or disassembled components for any remaining soil or disinfectant solution. If necessary, repeat the cleaning and disinfection. – Visually inspect the device or disassembled components for damage. In case of any damage to one or more parts, replace them with new ones.
			x	– Put spring back in its initial position. – Then push the release button carefully back into its initial position (see picture). – Verify that the release button springs back in position if pushed down to verify correct placement. – Re-insert silicone inlays in the handle. – Insert the axles pieces into Thopaz ⁺ back cover. – Slide handle into recess in axles and carefully tighten the screws using a Philips screw driver (see picture). – Verify that the bed holder is mounted securely and can be opened and closed.
				– To reassemble the docking station, re-install the spring on the shaft of the push button and re-insert the push button in slotted hole in the docking station. Next, hold the push button down and re-install the clamp and finally tighten the screw.
Storage	x		– Always store device in a dry, clean, and dust free environment.	
Transport	x	x	x	– Reprocess the device before sending it in for service. If this is not possible or can only be done in parts, the package shall indicate the potential biohazard. Local procedures and guidelines apply.

14 Warranty, maintenance and checks

Warranty

Medela AG warrants the device will be free from defects in materials and workmanship for a period of 2 years from the date of delivery ex works. Faulty material will be replaced free of charge during this period if not resulting from abuse or misapplication. This will not apply to parts subject to wear and tear in use. To ensure compliance with this warranty as well as optimum service from Medela products, we recommend the exclusive use of Medela accessories with our appliances.

In no event shall Medela AG be liable for claims which exceed the scope of warranty described including liability for consequential damages, etc. The right to the replacement of faulty parts will not be recognized by Medela if any work has been carried out on Thopaz+ by unauthorized persons. This warranty is subject to the appliance being returned to a Medela service centre.

Maintenance

Thopaz+ is a maintenance free device. Each time the device is switched on, a self-test is carried out which checks the internal functions of Thopaz+. A "beep" indicates that Thopaz+ successfully passed the self-test and that the device is switched on. Safe operation is ensured by the functional check including canister and tubing to be carried out before each use. The safety instructions must be observed.

Routine check

Medela recommends subjecting Thopaz+ to a routine check once a year. The work to be carried out and a protocol template in accordance to EN/IEC 62353 is given in Appendix B in these instructions. No additional calibration is required as long as Thopaz+ is used according to the instructions for use.


Instructions for FPT test

**Important: Patient must not be connected!
Thopaz+ must be disconnected from the mains
(battery operation) and turned off.**

Vacuum gauge
REF 077.1456


REF 079.0021

REF 079.0016

(Tube e.g. \varnothing 12x7 mm, not included in the scope of delivery)

1. Connect the tubing [REF 079.0021](#) and the 0.8l canister [REF 079.0016](#) to Thopaz+. For guidance to connect the system refer to chapter 7 – Preparation for use.
2. Connect vacuum gauge with tube to the patient connector.

3. 
 - 3.1 Press and hold [▶] and [■] simultaneously. Press additionally [⏻].
 - 3.2 The safety test starts.
 - 3.3 Follow the instructions on the display.

Press simultaneously

Safety check

There are no prescribed or recommended safety checks to be carried out.

Rationales

- The engineering of the Medela Thopaz+ has been verified by independent test institutes for compliance with the IEC 60601-1 standard. Copies of such certificates may be obtained on request. Because of the method of construction used, Medela does not expect electrical safety to be affected at any time during the life of the product – provided that Thopaz+ is repaired throughout its service life strictly and exclusively by Medela authorized service centers and that it is used properly in accordance with the intended use.
- For Thopaz+ as a device in protection class II (IEC 60601-1), there is no protective earth conductor; protective earth related tests (protective earth resistance, earth leakage current etc.) are not applicable.
- Thopaz+ enclosure is made entirely of insulating material. Tests of the touch current using common measuring instruments will therefore not reveal measurable values.
- The applied parts of Thopaz+ are connected to Thopaz+ via non-conductive vacuum hoses, canisters and tubings. Even when suctioning a conductive fluid (9g/L NaCl; IEC 60601-1) until the overflow protection device activates, measurements of the patient leakage current using common measuring instruments will not reveal measurable values.
- Thopaz+ systems do not have patient connections and functional earth conductors; tests related to those connections and currents are not applicable.

15 Disposal

Separate the parts of Thopaz⁺ system and dispose of according to the following instructions.

Sterile canister, tubing and caps

Please take care that you dispose of Thopaz⁺ sterile accessories in accordance with the hospital disposal guidelines.

Non-sterile general accessories

The parts are made of plastics that are not harmful to the environment when disposed of as household waste. Recycle or dispose of according to local regulations.

Pump and electrical parts

- The suction pump has to be disposed of in accordance with the European directive 2012/19/EU WEEE.
- Do not dispose of electric or electronic equipment together with unsorted municipal waste, rather collect it separately.
- In the European Union/ Switzerland/ UK the manufacturer or its vendor must take back waste equipment. Other countries may have similar collection and recycling systems. Please respect the relevant state laws and rules in your country for the disposal of electrical and electronic equipment.
- The separate collection and recycling of your waste equipment at the time of disposal will help conserve natural resources and ensure that it is recycled in a manner that protects human health and the environment.

Lithium-ion batteries

- The motor unit contains lithium-ion batteries which bear risk of fire, explosion and burns, if disposed of improperly.
- Do not disassemble, damage, short-circuit, heat above 80 °C (176 °F), incinerate or dispose of in fire.
- Dispose of in accordance with European Directive 2012/19/EU WEEE and local requirements.

16 Accessories

WARNINGS

Thopaz⁺ was validated in combination with the accessories listed in Appendix A. For a correct and safe operation use Thopaz⁺ with these accessories only. Further information is supplied with the individual accessory.

17 Technical specifications

low vacuum,
-10 kPa/-75 mmHg/
-102 cmHzO/-100 mbar
Tolerance: ± 15%

Measured at 0 m, atmospheric pressure:
1013.25 hPa
Please note: vacuum levels may vary depending on location (meters above sea level, atmospheric pressure and temperature).

≤ 5 l/min

1.1 kg/2.4 lbs

Power supply
Model: MSA-C2500IS12.0-30C-ZZ
IEC: 60601-1
Input: 100-240 VAC, max 0.8 A,
50/60 Hz
Output: 12 VDC, 2.5 A


Pump
12 VDC, 20W


ISO 13485
CE (2017/745/EU), IIb

HxWxD without canister
228 x 172 x 97 mm

Transport/Storage Conditions

Operating Conditions

IP33

Battery model

Rechargeable Lithium-ion battery

7.2 VDC 2.42 Ah
Weight 100 g ± 5g

HY-LINE A
Hochstrasse 355
8200 Schaffhausen
Switzerland
www.hy-line-group.ch
info@hy-line-group.ch
Made in Taiwan

18 Meaning of symbols

The following symbols may be present on the device labeling or in information accompanying the device. Refer to the device labeling or accompanying information for the applicable symbols.

	REF	MD	SN
General safety alert symbol, points to information related to safety.	Indicates the part number of the device.	Indicates the item is a medical device.	Indicates the serial number of the device.

			
Defines a temperature range (e.g. for operation, transport or storage).	Indicates the range of atmospheric pressure to which the medical device can be safely exposed.	Defines a relative humidity range (e.g. for operation, transport or storage).	Read and follow the instructions for use.

R_x only			
Indicates a Prescription Device U.S. Federal law restricts this device to sale by or on the order of a physician (for US only).	Identifies the manufacturer.	Contains fragile goods. Handle with care.	Keep away from rain. Keep in dry conditions.

			
Keep away from sunlight.	Indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.	MRI Unsafe keep away from magnetic resonance imaging (MRI) equipment.	Do not dispose of electric/ electronic devices together with unsorted municipal waste (dispose of the device in accordance with local regulations).

			
Indicates the location of the mute button.	Indicates a Class II electrical equipment.	Indicates alternating current.	Indicates the date of manufacturing.

LOT			
Indicates manufacturer's batch code.	Indicates the quantity (x) of individual devices in pack.	Indicates that the device should not be used after the date shown.	Indicates a single use device. Do not reuse the device.

STERILE EO	CE 0123	§	
Indicates the device is sterilized using ethylene oxide.	Indicates compliance with the EU requirements concerning medical devices.	Indicates the legal specifications of the system.	Indicates to consult instructions for use.

		IP33	
Indicates safety related tip.	Indicates the class of the system related to electrical safety.	Indicates the protection against ingress of solid foreign objects and against harmful effects due to the ingress of water.	Indicates a type CF applied part.

			
Indicates that interferences may occur in the vicinity of equipment marked with this symbol.	Indicates do not use the device if package is damaged.	Indicates that the package is capable of being recycled.	Indicates the maximum vacuum level of the system.
			
Indicates the flow levels of the system.	Indicates the electrical specifications of the system.	Indicates direct current.	Indicates the weight of the system.
	pcs AC		
Indicates the dimensions (h x w x d) of the system.	Indicates number of items.	Indicates air leak.	Indicates fluid.
			
Indicates Patient number.	Indicates Thopaz+ is out of axis and fluid cannot be measured.	Indicates to press the button more than 3 seconds.	Indicates that battery is empty.
			
Indicates that battery is almost empty.	Indicates that battery is being charged (moving bars).	Indicates that battery is fully charged and connected to the mains.	Indicates fluid cannot be measured.
			
Indicates not enough fluid for valid fluid measurement.	Indicates the country of manufacture of the product.	Indicates recycling of the packaging material defined with the code "XX" and the abbreviation "YYY".	Indicates a single sterile barrier system.
	UDI		#
Indicates a single sterile barrier system with protective packaging outside.	Indicates a carrier that contains unique device identifier information.	Indicates this side up.	Indicates model number.
		EC REP	EU REP
Indicates the entity importing the medical device into the locale.	Indicates the entity distributing the medical device into the locale.	Indicates the authorized representative in the European Union.	Indicates the authorized representative in the European Union.

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3.

The diagram shows the control panel of the Thoraz+ device. It features a central screen with the text 'Thoraz+' and 'Cardiorespiratory Therapy'. To the right of the screen are several control buttons: a power button (indicated by a '1' and a power symbol), a play button (indicated by a '2' and a right-pointing triangle), a left-pointing triangle, a right-pointing triangle, and a stop button (indicated by a '2' and a square with a diagonal line). Below the screen, there are two arrows pointing outwards, also labeled with '1'. A bracket at the bottom indicates that the power button and the two arrows should be pressed simultaneously.

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19 Technical documentation

WARNINGS

Do not use other accessories than those specified or sold by the manufacturer as replacement parts for internal components as it may result in increased emissions or decreased immunity of the Thopaz+ device.

EMC

Thopaz+ is EMC-tested in conformity with the requirements of IEC 60601-1-2:2014/AMD1:2020 Edition 4.1 according clause 7 and 8.9. Thopaz+ is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the EMC information. Thopaz+ does not have an essential performance as defined in IEC 60601-1.

Electromagnetic compatibility

Electromagnetic emissions

Thopaz+ is only approved for the following electromagnetic environments: professional healthcare facility environment and home healthcare environment.

Emision Tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	Thopaz+ uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference with nearby electronic equipment.
RF emissions CISPR 11	Class B	Thopaz+ is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

WARNINGS

- Thopaz+ should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, Thopaz+ should be observed to verify normal operation in the configuration in which it will be used.
- Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies, RFID can affect the Thopaz+ pump and should be kept at a distance at least 30 cm away from the device.
- Use of accessories or cables other than those provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Electromagnetic immunity

Thopaz+ is only approved for the following electromagnetic environments: professional healthcare facility environment and home healthcare environment.

Immunity Tests	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV line-to-earth	± 1 kV differential mode ± 2 kV line-to-earth not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<0% U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T for 1 cycle 70% U_T for 25 cycles at 50 Hz single phase: at 0° and for 30 cycles at 60 Hz single phase: at 0° <0% U_T for 250 cycles at 50 Hz and for 300 cycles at 60 Hz	<0% U_T for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% U_T for 1 cycle 70% U_T for 25 cycles at 50 Hz single phase: at 0° and for 30 cycles at 60 Hz single phase: at 0° <0% U_T for 250 cycles at 50 Hz and for 300 cycles at 60 Hz	Mains power quality should be that of a typical commercial or hospital environment. If the user of Thopaz+ requires continued operation during power mains interruptions, it is recommended that Thopaz+ be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels of a typical commercial or hospital environment.

Proximity magnetic fields IEC 61000-4-39	8 A/m 30 kHz – CW	8 A/m 30 kHz – CW	Magnetic field intensity should be that of a typical or commercial or hospital environment.
	65 A/m 134.2 kHz – PM 2.1 kHz	65 A/m 134.2 kHz – PM 2.1 kHz	
	7.5 A/m 13.56 MHz – PM 50 kHz	7.5 A/m 13.56 MHz – PM 50 kHz	
<p>NOTE U_T is the a.c. mains voltage prior to application of the test level. CW: Continuous Wave PM: Pulse Modulation</p>			

Electromagnetic immunity

Thopaz+ is only approved for the following electromagnetic environments: professional healthcare facility environment and home healthcare environment.

Immunity Tests	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3Vrms 0.15–80MHz 6Vrms in ISM and amateur radio bands between 0.15 and 80MHz	3Vrms 6Vrms	<p>Portable and mobile RF communications equipment should be used no closer to any part of Thopaz+, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance</p> $d = 1.2 \sqrt{P}$ $d = 0.35 \sqrt{P}$ <p>80MHz to 800MHz</p> $d = 0.7 \sqrt{P}$ <p>800MHz to 2.7GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	10V/m 80MHz to 2.7GHz	10V/m	

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

NOTE 3 Proximity fields from RF wireless communications equipment were tested according to Table 9 of IEC 60601-1-2:2014/AMD1:2020

^a Field strengths from fixed RF transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which Thopaz+ are used exceeds the applicable RF compliance level above, Thopaz+ should be observed to verify normal operation. If abnormal operation is observed, additional measures may be necessary, such as reorienting or relocating Thopaz+.

^b Over the frequency range 150kHz to 80MHz, field strengths should be less than 10V/m.

Table of frequencies

Table of frequencies of portable and mobile transmitters for which the recommended separation distance is 30 cm (12 inches):

Band (MHz)	Service
380–390	TETRA 400
430–470	GMRS 460, FRS 460
704–787	LTE Band 13, 17
800–960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5
1 700–1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS
2 400–2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7
5 100–5 800	WLAN 802.11 a/n

WARNING

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm to any part of the Thopaz+ pump including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Appendix A

REF

- 079.1000 Thopaz+ digital cardiothoracic drainage system with Instruction for use 200.6841
- 101046671* 079.1002 Thopaz+ digital cardiothoracic drainage system with Instruction for use 200.6842
- 079.1003 Thopaz+ digital cardiothoracic drainage system with Instruction for use 200.6843/101034256
*not CE marked

Disposables

- 079.0011 Thopaz canister 0.3l, antifoaming*
- 079.0012 Thopaz canister 0.3l with solidifier, antifoaming*

- 079.0016 Thopaz canister 0.8l, antifoaming*
- 079.0017 Thopaz canister 0.8l with solidifier, antifoaming*

- 079.0018 Thopaz canister 2l, antifoaming*
- 079.0019 Thopaz canister 2l with solidifier, antifoaming*

- 079.0021 Thopaz tubing single*
- 079.0022 Thopaz tubing double*
- 079.0023 Thopaz tubing single, small connector*
- 079.0024 Thopaz tubing double, small connector*
- 079.0025 Thopaz tubing single, large connector*
- 079.0026 Thopaz tubing double, large connector*

- 101039326 Thopaz tubing extension, 1 line*

- 101039327 Thopaz tubing extension, 3 lines*

Accessories

079.1006 Carrying strap Thopaz+

079.0037 Docking station

079.0038 Adapter docking station for Thopaz+ canister 2l

079.0034 USB cable

079.0039 Thopaz sealing cap*

079.0036 Holder with standard rail

079.0040 Universal holder with standard rail

077.1456 Vacuum gauge

Documents

200.6841 Thopaz+ Instr. for use (EN/DE/FR/IT/NL/ES/PT/HU/PL)
200.6842 Thopaz+ Instr. for use (EN/DA/FI/NO/SV/CS/HR/SR/EL)
200.6843 Thopaz+ Instr. for use (EN/JA/KO/RU/TR/ET/LT/LV/SK)
101034256 Thopaz+ Instr. for use (EN/ZH)

Spare parts

087.0059 Thopaz+/Liberty wall mains adapter/charger
077.0117 Screw covers
077.0145 Feet
099.0343 Seal orange
077.0124 Cover for ports
099.0361 Standard rail holder
099.0348 Release button

For additional spare parts, please refer to the service manual.

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Medela Medical AB
Medela Schweiz AG
Medela UK Ltd.
Medela LLC