



Double Electric Breast Pumps	Pump In Style® Pro Pump In Style® Pro+	Pump In Style® Pro+ Hands-free	Swing Maxi™ Hands-free	Freestyle™ Hands-free
Description	<ul style="list-style-type: none"> ✓ Similar pumping pattern to Symphony®'s clinically proven Maintain mode ✓ Traditional bottles and shields ✓ Rechargeable battery (Pro+ only) ✓ Delivering 11.8% more milk* ✓ 21 and 24mm PersonalFit™ PLUS breast shields 	<ul style="list-style-type: none"> ✓ Rechargeable battery ✓ Similar pumping pattern to Symphony®'s clinically proven Maintain mode ✓ Hands-free mobility ✓ Delivering 11.8% more milk* ✓ 21 & 24mm in bra collection cups 	<ul style="list-style-type: none"> ✓ Mobile pumping ✓ Rechargeable battery ✓ Connected ✓ Lanyard for extra mobility ✓ Delivering 11.8% more milk* ✓ 21 & 24mm in bra collection cups 	<ul style="list-style-type: none"> ✓ Compact ✓ Rechargeable battery ✓ Connected ✓ Lanyard for extra mobility ✓ Delivering 11.8% more milk* ✓ 21 & 24mm in bra collection cups
Technology	2-Phase Expression® technology Added massage technology for your comfort and more milk	2-Phase Expression® technology Added massage technology for your comfort and more milk	2-Phase Expression® technology	2-Phase Expression® technology
mmHg*(vacuum) cycles per minute	-50 to -295 mmHg ^{1,3} 35-120 cpm	-50 to -295 mmHg ^{1,3} 35 to 120 cpm	-45 to -245 mmHg ^{1,2} 45-111 cpm	-45 to -270 mmHg ^{1,2} 45 – 111 cpm
Number of levels	16 ●●●●●●●●●●●●●●	16 ●●●●●●●●●●●●●●	9 ●●●●●●●●●	9 ●●●●●●●●●
Display / controls	16 level indicators 4 easy to use push buttons Pumping mode indicator for stimulation vs expression mode Level and battery life (Pro+ only) indicators	16 level indicators 4 easy to use push buttons Pumping mode indicator for stimulation vs expression mode Level and battery life indicators	4 easy to push buttons Display featuring in-app visibility to session timer, battery status, vacuum level, Bluetooth® connection, etc.	Touch screen display Illuminated: session timer, battery status, vacuum level, Bluetooth® connection, etc.
Main power	Plug-in adaptor global compatibility 100 – 240 V	Plug-in adaptor global compatibility 100 – 240 V	USB chargeable Plug-in adaptor global compatibility 100 – 240V	USB chargeable global compatibility 100 – 240 V
Battery options	Pro: Wall powered Pro+: Integrated rechargeable	Integrated rechargeable	Integrated rechargeable	Integrated rechargeable
Noise level				
Connectivity	–	–		
Breast Shields Sizes <small>Additional sizes available, (sold separately)</small>	 PersonalFit™ PLUS shields Oval shape for better fit and 105° angle for greater comfort and milk flow	 In-bra collection cups Anatomic, discreet design and 105° angle for greater comfort and milk flow	 In-bra collection cups Anatomic, discreet design and 105° angle for greater comfort and milk flow	 In-bra collection cups Anatomic, discreet design and 105° angle for greater comfort and milk flow
Breast Milk Collection Containers	(2) 5oz, 150mL bottles	(2) 5oz, 150mL Collection Cups	(2) 5oz, 150mL Collection Cups	(2) 5oz, 150mL Collection Cups

Content appearance may vary from pictures. Contents included may vary by breast pump version.

*Clinical study, data on file.
 References 1 Maximum Vacuum in Expression (when single or double pumping). 2 Under the typical conditions, the Swing Maxi and Freestyle breast pumps are capable of providing vacuum levels from -45 to -245 mmHg. Under the foreseeable limits of use conditions which produce peak vacuum of the greatest magnitude, the pump could produce vacuums of -270 mmHg. 3 Under the typical conditions, the Pump In Style breast pump is capable of providing vacuum levels from -50 to -240 mmHg. Under the foreseeable limits of use conditions which produce peak vacuum of the greatest magnitude, the pump could produce vacuums of -295 mmHg.