

Supplemental Nursing System (SNS)

The Supplemental Nursing System is a nursing device specifically designed to help mothers and parents to breastfeed / chestfeed by supplying additional milk to the baby whilst they are being fed at the breast / chest.

- Facilitates bonding and skin-to-skin contact
- Promotes milk supply by stimulating the breast through sustained sucking
- Supports adoptive, surrogate, same sex, transgender and non-binary parents to create a breastfeeding / chestfeeding experience
- Supports baby's oral development of sucking skills to create vacuum
- Provides milk supplementation without the use of bottles and teats

Optimised retail packaging with hook



What's included



1x Container & Lid



1x TwistLok Lid



1x TwistLok Tube



1x Ring



1x Cleaning Aid



1x Clamp



1x Clip

1x Instruction for use



Advanced system

- Anti-leak TwistLok Lid prevents leakage and milk loss
- Advanced venting system ensures consistent milk flow and comfortable feeding
- Silicone tubing with soft, rounded edges - gentle on baby's mouth



Easy to use

- Versatile container and lid, ideal for milk preparation, storage, transport and feeding
- Clip securely fastens device to nursing bra, top or pillow
- Simple "on/off" flow control



Easy to clean

- Dedicated cleaning aid to help effectively clean the tube over repeated daily use



Reusable



Swiss design

Product Name	Item number	Languages	Retail box (1 unit)			Shipper			
			EAN/ UPC code	Dimensions (LxWxH in mm)	Weight (kg)	EAN/ UPC code	Qty/Case	Dimensions (LxWxH in mm)	Weight (kg)
Supplemental Nursing System	101043648	EN ES FR	020451436487	117 x 62.5 x 131		20020451436481	6	271 x 195 x 130	
Supplemental Nursing System	101043649	DE, FR, IY, NL, ES, PT, PL, CS, SK, HR, SR, SL	7612367080439			7612367080446			
Supplemental Nursing System	101043650	EN, DA, SV, NO, EL, ET, RU, HU, RO, JA, ID (IFU only)	7612367080453			7612367080460			

To order or for more information, contact Medela Customer Service at customer.service@medela.com

Scan the QR code for further SNS information

